
Chapter 1

OMA-206-P Introduction

Introduction:

The OMA-206-P is a portable diode-array fiber-optics spectrophotometer, continuously measuring the absorbance spectra of the given sample. It is constructed of a light source, a dispersion device for selecting a particular wavelength, a sample area, and a detector. The OMA uses a xenon light source, a holographic grating to disperse the light according to its wavelength, a cuvette that contains the sample to be measured or the blank (zero) sample, and a diode array detector which converts the light intensity into electrical signal.

A full spectrum is continuously measured and analyzed, via a multi wavelength method, to give the composition of the stream.

The OMA-206-P is designed to address on-line applications that require cost effective, accurate, and continuous concentration monitoring. A simple touch-screen based user interface is ideal for single and multi component analysis and allows for quick setup and calibration.

Wavelength ranges from 200-800 nm are continuously and instantaneously monitored by a diode array spectrophotometer with 1024 elements. A nonvolatile memory section of the analyzer contains the calculation method, wavelength ranges, and the name of the analyte to be displayed on the screen.

The OMA-206-P is housed in a weather resistant ruggedized case. The case is IP67, MIL C-4150J, STANAG 4280, and ATA 300 certified. Two fiber optic cables are used to conduct light to and from the cuvette, where the samples are measured.

The zeroing process is manual. The user blanks the analyzer on a cuvette full of zeroing solution then inserts a sample. A button on the touch screen is pushed to store the zero sample and a the analysis.

The measurements and analysis are continuous.

Note: Instructions in this manual regarding Auto zero should be ignored since the OMA-206-P does not include a sampling system.



OMA-206-P Specifications

Analyzer

Analyzer type

OMA-206-P portable spectrophotometer

Wavelength range

190-800nm

Slit Width

1nm

EP resolution

>1.6 toluene in hexane,
ratio absorbance at 269nm / 266nm

Wavelength Accuracy

<+-0.5nm (NIST 2034)

Wavelength reproducibility

<+-0.04nm (NIST 2034)

Photometric noise

<0.002AU 32 scans at 0 AU at 250nm

Photometric stability

<0.002AU/h at 0AU at 340nm +-10C

Operating conditions

Temperature

0⁰ to 55⁰ C (32⁰ to 130⁰ F)

Display

1/4" VGA , NEMA 4 touch screen

LCD display (340x240 pixels)

Power

80 to 240 Volts AC; 40 to 60 Hz ; 20 W


Analyzer Dimensions

16"(H)x14"(W)x8"(D)(41x36x20cm)


Area Classifications

IP67, MIL C-4150J, STANAG 4280, and ATA 300 certified

REVISIONS			
REV	DESCRIPTION	DATE	APPROVED
0	ELECTRICAL DIAGRAM		A.R.



OMA-206-P PANEL

<div><div><div></div><div></div><div></div></div><div><div>Applied Analytics, Inc.</div><div>4 Clock Tower Place, Suite 420 Maynard, MA 01754</div></div></div>		TITLE		ELECTRICAL DIAGRAM, OMA-206-P PORTABLE ANALYZER		SIZE A	DWG. NO. ED-334-0002	REV. 0
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE: DECIMAL ANGLES 63° XX.015 XXX.01 XXX.005		COMPANY CONFIDENTIAL THE DOCUMENT AND THE DATA DISCLOSED HEREIN OR HEREWITH IS NOT TO BE REPRODUCED, USED, OR DISCLOSED IN WHOLE OR IN PART TO ANYONE WITHOUT THE WRITTEN PERMISSION OF APPLIED ANALYTICS, INC.		APPROVALS		DATE	SCALE XX	
THIRD ANGLE PROJECTION		DO NOT SCALE DRAWING		DRAWN E.M.			SHEET 1 OF 1	
		BREAK ALL SHARP EDGES		CHECKED V.J.				
REMOVE ALL BURRS AND BREAK SHARP EDGES		MATERIAL		APPROVED A.P.				
INDIVIDUALLY PACKAGE AND LABEL SMALL PARTS WITH THE DRAWING NUMBER								