Environma Water Vapour Permeation Analysers Stion Analysing

nemation

ency Ultraviolet Gas Production owdered Metals d Environments

s Manufacturing Carbon Refining

Food Packaging s ■ Glove Boxes on Beam ■ R & D ■ Fermentation d Environments Manufacturing ssel Blanketing

tion Analysing ncy Ultraviolet Gas Production Ndered Metals

Environments anufacturing el Blanketing

S Production ered Metals

ufacturing

ve Boxes

ıfacturing

Analysing Ultraviolet

; Production



Modular systems for precision water vapour analysis of packaging film barriers



Applications

Barrier films PET bottles Containers Closures Flexible Pouches Bags

Features & Benefits

- **Analytical Systems Manufactured** traceable to NIST.
- System validation with certified gas or film for speed and convenience.
- Over 25 yrs experience of Proprietary Coulometric P₂O₅ sensor.
- Absolute moisture measurement No calibration required.

- Widest measuring range in the market providing research grade flexibility.
- Electronic flow, temperature and humidity control for ultimate responsiveness and repeatability.
- Intuitive Windows based software.
- No archaic liquid coolants, catalysts or special gas mixtures required.

Conforms to: **ASTM F-1249*** ISO 15106-3 ISO 15105-2 DIN 53122-2

Environments lanufacturing sel Blanketing tion Analysing ncy Ultraviolet Gas Production owdered Metals d Environments Manufacturing carbon Refining Food Packaging ■ Glove Boxes n Beam ■ R & D Fermentation l Environments Manufacturing isel Blanketing tion Analysing cy Ultraviolet as Production dered Metals nvironments anufacturing l Blanketing Production red Metals

Manufacture traceable to NIST

All Systech analysers are certified traceable to NIST. In addition, analytical performance is validated using NIST certified gases and NIST traceable films. A set of validation films and a spare P_2O_5 sensor comes as standard with all of the water vapour permeation analysers.

P₂O₅ Sensor Technology

The two most common sensor types for measuring moisture are IR (infra-red) and P_2O_5 (phosphorous pentoxide). The P_2O_5 sensor is more sensitive and stable than IR and does not require calibration. P_2O_5 is the primary method for absolute moisture measurement. Systech have over 20 years experience of using P_2O_5 sensor technology.

*The Systech 7000 series analysers comply with ASTM standard F-1249 with the exception of the sensor technology. The standard relates to an infra red sensor whilst the Systech analyser uses a coulometric sensor - a dedicated method of moisture analysis.

Widest measuring range

The 7000 series offer you the widest measurement range in the market providing research grade flexibility.

- Measurement range of 0.002-1000 g/m²/day with masking.
- Quality Assurance orientated speed and agility.
- Up to five expansion modules available to increase testing throughput.



Precision control

These analysers offer precision electronic, temperature and humidity flow control providing ultimate responsiveness and repeatability.

- Test gas and carrier flow gas controlled by premium electronic mass flow controllers.
- Widest sample temperature range available of 5 to 50°C.
- Accurate relative humidity range from 0% to 90%.

Laboratory Testing Services

Our test laboratory will perform your Permeation Testing and Headspace Analysis. Whether you are developing innovative materials and packages or validating that your supplier is meeting specification. We can exceed your expectations with:

- Competitive Prices
- Fast Turnaround

- Independent non-biased results
- 25 Years Experience

Software `


The intuitive Windows based software offers:

- Easy input and recall of operating parameters and test protocols.
- User-friendly data tracking, searches, sorts, storage and output capabilities.
- Graphical representation of measurement data in real time.
- Auto-stop feature stops test when samples have reached equilibrium or by user entered elapsed time value.
- Complete system diagnostics.


Start up



Set up



Actual data




Principle of Operation

Utilising our proprietary sensor technology to detect water vapour transmission rates, samples are clamped or attached to a diffusion chamber. Wet nitrogen is then introduced into the upper half of the chamber while a moisture-free carrier gas flows through the lower half.


Molecules of water diffusing through the sample into the lower chamber are conveyed to the sensor by the carrier gas.

Sample Test Chamber



This allows a direct measurement of the water vapour without using complex extrapolations. Water vapour transmission rate of the test sample is displayed as either g/m²/day or g/100in²/day.

P2O5 Sensor



To achieve an absolute measure, the technology draws upon a fundamental principle of physics.

The phosphorous pentoxide (P₂O₅) moisture sensor consists of a dual platinum winding formed around a quartz tube.

The change in the resistance across the windings creates a change in the measured current. According to Faraday's Law this is directly proportional to the amount of moisture in the gas stream.

ontact Leng inealing **=** V amics Comb Oxygen Defi s 🔳 High Puri ■ Alloys and lity ■ Control ı ■ Contact Le keting ■ Hyd on Analysing ass/Fibre Opt uction ■ Elect Ndered Metal ity ■ Control ■ Contact Ler Annealing ■ \ ramics Comb Oxygen Defici ■ High Purity Alloys and Pa ■ Controlled ontact Lens aling ■ Ves s Combust n Deficien and Pow ntrolled Er tics G als ■ Fe ■ Vesse ombustio

Deficienc Purity G

911012

Systech 7000 Series - Water Vapour Permeation Analysers

Systech's range meets the requirement for the testing of any application.



Systech 7002 Water Vapour Transmission Analyser delivers the same high performance as the 7001 but with an extended measurement range for more demanding applications.

Technical Specifications

Measurement Range

7001 Unmasked 0.002 to 10 g/m²/day 0.02 to 70 g/m²/day Masked

7002 Unmasked 0.002 to 70 g/m²/day Masked 0.02 to 1000 g/m²/day

Test Temperature Range 5°C to 50°C (41°F to 122°F)

Test RH Range 0 to 90% RH

Sample Size 50cm², adaptors available for smaller samples

Operating Conditions Standard laboratory environment

Power Requirements 110/220 VAC 50/60Hz

Supply Pressure 25 - 45 PSI regulated

Gas Fittings 1/8 in. Swagelok (supplied)

Enclosure Epoxy coated heavy gauge steel

Dimensions 533 x 533 x 305 (mm)

Weight 23.6kg

Options

Expansion Modules 7011 Available for simultaneous WVTR

analysis of up to 12 samples

Environmental Chambers For finished package testing

Systech Instruments have over 25 years experience of providing analysis solutions for a wide range of industries. From our manufacturing plant in the UK we produce gas analysers for industrial process industries, headspace analysers for monitoring gas flushing of food products, and our range of permeation analysers.



17 Thame Park Business Centre, Wenman Road, Thame, Oxfordshire OX9 3XA

Tel: +44 (0)1844 216838 Fax: +44 (0)1844 217220 E-mail: advice@systech.co.uk www.systechinstruments.com



Illinois Instruments, inc 2401 Hiller Ridge Road Johnsburg, Illinois 60050

Tel: +1 815 344 6212 Fax: +1 815 344 6332

E-mail: sales@illinoisinstruments.com www.illinoisinstruments.com

ontact Lens I inealing ■ Ve amics Combu Oxygen Defici s
High Purity ■ Alloys and P lity ■ Controlle ■ Contact Len keting ■ Hydro on Analysing **•** ass/Fibre Optio uction ■ Electro Ndered Metals ity ■ Controlle ■ Contact Lens Annealing ■ Ve ramics Combus Oxygen Deficie ■ High Purity Alloys and Pov ■ Controlled | ontact Lens M aling ■ Vesse s Combustic n Deficiency h Purity Gas and Powd ntrolled Env Lens Mani ■ Vessel I ombustion Deficiency

Purity Gas

pallonn